

Reference these catalogue numbers to order:

B150-048 - Bioburden Filtration Funnels, 150ml, 48 funnels per case

B250-048 - Bioburden Filtration Funnels, 250ml, 48 funnels per case


Part number	Name	Shelf Life	Storage conditions
SKU-B150-048	BB Filtration Kit, 150ml	18 Months	Room temperature
SKU-B250-048	BB Filtration Kit, 250ml	18 Months	Room temperature

**Product Description:** The filtration kit consists of a black 0.45-micron pore size nitrocellulose ester membrane held in a plastic base below a removable filtration funnel. Traditionally the 150mL filtration funnels are used for Bioburden testing and the 250mL for Water testing. The black filter used with our funnels is optimized for compatibility with the Growth Direct® System.

**Irradiation:** Our B150-048 and B250-048 are produced in an ISO9001-certified facility and sterilized in accordance with Rapid Micro Biosystems® Quality Systems. This is a Gamma irradiation treated product at 15-40kGy, packaged as 8 sets of 6 units in poly sleeves stored in a cardboard box.

**Material:** Polypropylene, Color: Natural, Material Grade: PRO-FAX, PF511.





RMB strives to provide the highest quality products and services to our customers. If you have any questions regarding this informational document, please contact Customer Service at [support@rapidmicrobio.com](mailto:support@rapidmicrobio.com). If interested in placing an order with Rapid Micro Biosystems®, please reach out to [orders@rapidmicrobio.com](mailto:orders@rapidmicrobio.com).