CERTIFICATE OF QUALITY

Catalog No: EMVL-100 Lot No: 106252 Expiry: 30JUN2022
Storage: Room temperature

PRODUCT DESCRIPTION:

The Cassette Cover consists of a transparent non-fluorescing plastic window that allows images to be taken by the Growth Direct of the membrane on the enclosed media surface.

MANUFACTURE:

The product was manufactured in an ISO certified facility. The sterilization treatment was performed in an ISO certified facility.

QUALITY RELEASE CRITERIA:

The quality control data given in this document represent the quality of the released lot as tested.

PHYSICAL TESTS:

1. Optical background	Pass
2. No fluorescent debris	Pass

This is sterilized as indicated on the product label, packaged and sealed in a Tyvek/poly pouch and stored in a cardboard box.

Released by:
Quality Assurance

Date: 21 Sep 2020

Rapidmicro biosystems
1001 Pawtucket Blvd
Lowell, MA 01854
978-349-3200

This product is not for In Vitro diagnostic use.
Country of origin- USA