CERTIFICATE OF QUALITY

Catalog No: EMVL-100 Lot No: 96140 Expiry: 31JAN22
Storage: Room temperature

PRODUCT DESCRIPTION:

The Cassette Cover consists of a transparent non-fluorescing plastic window that allows images to be taken by the Growth Direct of the membrane on the enclosed media surface.

MANUFACTURE:

The product was manufactured in an ISO certified facility. The sterilization treatment was performed in an ISO certified facility.

QUALITY RELEASE CRITERIA:

The quality control data given in this document represent the quality of the released lot as tested.

PHYSICAL TESTS:

1. Optical background	Pass
2. No fluorescent debris	Pass

This is sterilized as indicated on the product label, packaged and sealed in a Tyvek/poly pouch and stored in a cardboard box.

Released by: _____ Date: 07 JUL 2020
Quality Assurance

Rapidmicro biosystems
1001 Pawtucket Blvd
Lowell, MA 01854
978-349-3200

This product is not for In Vitro diagnostic use.
Country of origin- USA